



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sternberg Double Ended Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sternberg Double Ended Retractor

SKU RM-SDER-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
14 cm Baseline — U-shaped opposing wire ends	RM-SDER-1-1	\$24
37 cm Macro-Reach — U-shaped opposing wire ends	RM-SDER-1-2	\$34

PRODUCT DESCRIPTION

The Sternberg Double Ended Retractor is an essential oral, dental, and maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic displacement of the cheeks, lips, and oral mucosa during intraoral procedures. Designed with an open wireframe chassis, this instrument shields delicate labial commissures and buccal tissue planes from dental burs and surgical scalpels without crowding the operative field.

Built to provide maximum durability and backed by a lifetime warranty against rust and material defects, this Sternberg Double Ended Retractor series offers two operational structural metrics: a 14 cm Length wireframe construction baseline (\$24) for standard oral surgery and an extended 37 cm Length macro-reach configuration option layout (\$34) for deep pharyngeal and transoral access. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry configuration incorporates an asymmetrical double-ended design featuring U-shaped curved ends facing opposite directions with smooth rounded wire edges that prevent mucosal pinching.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Sternberg Double Ended Retractor to withstand repetitive hospital autoclave cycles, providing oral surgery centers, dental clinics, and facial plastic surgery suites worldwide with exceptional corrosion resistance, structural resilience, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)