



RETRACTOR MAKER

PRODUCT BROCHURE OF

Yancoskie Coronal Lift Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Yancoskie Coronal Lift Retractor

SKU RM-YCLR-9IN-1

\$119 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Yancoskie Coronal Lift Retractor	9" (22.9 cm) overall length, 150 mm wide rake, 6 counterbalance sharp prongs	\$119
Protective Sleeve / Box	Individual transit packaging	\$119

PRODUCT DESCRIPTION

The Yancoskie Coronal Lift Retractor is a specialized facial plastic and aesthetic surgical retractor engineered specifically to provide safe, expansive scalp flap elevation during coronal brow lifts, subgaleal forehead rejuvenation, and fronto-orbital reconstructions. Designed with an extra-wide rake head, this instrument allows surgeons to apply uniform manual counter-traction across the calvarium, elevating the scalp flap cleanly while safeguarding delicate supraorbital and supratrochlear neurovascular bundles from crush trauma.

Built to withstand high-frequency autoclave sterilization cycles without structural degradation, this Yancoskie Coronal Lift Retractor configuration features an overall tool length of 9 inches (22.9 cm, \$119). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure supreme resistance to environmental corrosion, the active working width incorporates a 150 mm wide prong interface. The prong architecture features a single-ended layout integrating 6 counterbalance sharp prongs specifically angled to grip dense galeal tissue securely without slipping or puncturing through the cutaneous surface. Its elongated handle assembly provides operating assistants with fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Yancoskie Coronal Lift Retractor to uphold strict clinical standards under high-intensity surgical illumination, providing aesthetic plastic surgery clinics, facial reconstructive suites, and craniofacial theaters worldwide with exceptional durability, corrosion resistance, and complete autoclave dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)