



RETRACTOR MAKER

PRODUCT BROCHURE OF

Carpal Tunnel Self-Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Carpal Tunnel Self-Retaining Retractor

SKU RM-CTSRR-1

\$144 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Carpal Tunnel Self-Retaining Retractor	4 ¼" (10.79 cm) length, blunt prongs, ratchet lock	\$144
Protective Sleeve / Box	Individual transit packaging	\$144

PRODUCT DESCRIPTION

The Carpal Tunnel Self-Retaining Retractor is a precision-engineered hand and microsurgical retractor designed specifically to provide hands-free, atraumatic wound margin retraction during carpal tunnel release surgery. By securing palmar cutaneous flaps and subcutaneous fat under steady, locked mechanical tension, this instrument establishes an unobstructed visual field over the flexor retinaculum while eliminating manual assistant fatigue.

Built specifically for the confined anatomy of the hand and wrist, this Carpal Tunnel Self-Retaining Retractor configuration features a compact 4 ¼-inch (10.79 cm) overall length (\$144). Its working ends are equipped with smooth, blunt-tipped prongs engineered to safely displace superficial soft tissues without risk of puncturing the median nerve, its recurrent motor branch, or superficial palmar vascular arches. The integrated mechanical ratchet lock prevents accidental movement during delicate ligament division.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Carpal Tunnel Self-Retaining Retractor to strict clinical tolerances, ensuring smooth ratchet action, ergonomic finger-ring leverage, and complete autoclave compatibility for plastic, hand, and orthopedic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)