



RETRACTOR MAKER

PRODUCT BROCHURE OF

Finochietto (Modif) Rib Spreader

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Finochietto (Modif) Rib Spreader

SKU RM-FMR-299-1

\$299 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Streamlined Horizontal Track Frame	Flat slide rail supporting deep access manual widening (120 mm spread)
2 (Fixed Pair)	Fixed Curved Blades	14 mm × 15 mm fenestrated lateral blades contoured to rib surfaces
1	Rotary Screw Lock Assembly	Tightening gear to secure optimum surgical retraction
2	Ergonomic Forceps Extension Arms	Lightweight tissue-spreading handles to decrease user fatigue
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Finochietto (Modif) Rib Spreader is a specialized thoracic surgical retractor engineered specifically to provide safe, gradual rib separation and chest wall exposure during anterior and lateral thoracotomy procedures. Designed with fixed fenestrated lateral blades, this modified instrument allows surgeons to apply calibrated mechanical spreading forces across delicate rib structures, preventing bone fractures while maintaining a stable, unhindered visual corridor into the pleural cavity.

Built to clinical framework specifications, this Finochietto (Modif) Rib Spreader configuration (\$299) features a streamlined horizontal track frame with a flat slide rail that supports deep access manual widening up to a 120 mm maximum expansion capability. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade dimensions incorporate 14 mm × 15 mm fixed curved blades contoured to rib surfaces, making them ideal for delicate bones and pediatric or slender adult thoracic cages. The rotary screw lock assembly tightening gear secures optimum surgical retraction, while two ergonomic forceps extension

arms decrease user fatigue.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Finochietto (Modif) Rib Spreader to withstand repetitive hospital sterilization cycles, providing cardiothoracic suites, trauma centers, and pediatric thoracic operating rooms worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)