



RETRACTOR MAKER

PRODUCT BROCHURE OF

Obwegeser Soft Tissue Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Obwegeser Soft Tissue Retractor

SKU RM-OSTR-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Obwegeser Soft Tissue Retractor	8 ¾" (22.2 cm) overall length	\$39
Protective Sleeve / Box	Individual transit packaging	\$39

PRODUCT DESCRIPTION

The Obwegeser Soft Tissue Retractor is a specialized maxillofacial surgical retractor engineered to provide clear, reliable operative exposure while protecting sensitive soft-tissue structures during oral and orthognathic procedures. Designed with a curved distal channel lip, this instrument securely holds back the cheek mucosa, lips, and mucoperiosteal flaps away from oscillating saws and surgical drills.

Engineered for precision cranio-maxillofacial interventions, this standard Obwegeser Soft Tissue Retractor configuration features an overall length of 8 ¾ inches (22.2 cm). Its anatomically contoured blade features upward or downward curved channel lips that engage bone borders smoothly, preventing slippage and shielding facial neurovascular structures during complex ramus and maxillary osteotomies.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish, Retractor Maker produces the Obwegeser Soft Tissue Retractor to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and plastic surgeons with an ergonomic, durable instrument for reconstructive facial surgery.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)