



RETRACTOR MAKER

PRODUCT BROCHURE OF

Middledorpf Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Middledorpf Retractor

SKU RM-MDR-28-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
7.0 mm × 14.0 mm	RM-MDR-28-1-1	\$49
22.0 mm × 20.0 mm	RM-MDR-28-1-2	\$54
30.0 mm × 26.0 mm	RM-MDR-28-1-3	\$59

PRODUCT DESCRIPTION

The Middledorpf Retractor is an internationally acclaimed surgical instrument engineered specifically to provide safe, non-traumatic soft-tissue displacement and muscular retraction during deep general, pelvic, and orthopedic procedures. Designed with an open fenestrated curved blade, this retractor allows surgical teams to displace muscle masses and fascial layers cleanly while allowing light and microscopic vision to pass directly through the blade window.

Built for complete procedural adaptability, this Middledorpf Retractor maintains a uniform total length frame of 28 cm that provides excellent reaching clearance in deep cavities. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers two handle structural styles: a classic Finger Ring Handle for positive leverage, or a lightweight Hollow Handle engineered to reduce manual strain during extended cases. Surgeons can select from three fenestrated blade dimension profiles: 7.0 mm × 14.0 mm (\$49), 22.0 mm × 20.0 mm (\$54), and 30.0 mm × 26.0 mm (\$59).

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Middledorpf Retractor to withstand repetitive hospital autoclave cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)