



RETRACTOR MAKER

PRODUCT BROCHURE OF

Flared Cobra Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Flared Cobra Retractors

SKU RM-FCR-15-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Left-Facing — Single Prong, 15" (38 cm)	RM-FCR-15-1-1	\$69
Left-Facing — Double Prong, 15" (38 cm)	RM-FCR-15-1-2	\$69
Right-Facing — Single Prong, 15" (38 cm)	RM-FCR-15-1-3	\$69
Right-Facing — Double Prong, 15" (38 cm)	RM-FCR-15-1-4	\$69

PRODUCT DESCRIPTION

The Flared Cobra Retractors are specialized orthopedic surgical instruments engineered specifically to provide powerful bone leverage and expansive soft-tissue retraction during primary and revision total hip arthroplasty. Designed with an expansive flared blade head that arches smoothly backward, these retractors allow surgical teams to displace large gluteal and thigh muscle masses away from the joint corridor without causing muscle contusions.

Built for deep pelvic and joint clearance, this Flared Cobra Retractors series features a unified overall tool length of 15 inches (38 cm, \$69) that keeps the assistant's hands safely outside the surgical corridor. Surgeons can customize the instrument across directional variations—Left-facing orientation or Right-facing orientation—paired with selectable tip configuration profiles including Single Prong for focused bone purchase or Double Prong for balanced dual-point stability along the acetabular rim.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures Flared Cobra Retractors to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing joint reconstruction and orthopedic trauma teams worldwide with unyielding clinical durability and performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)