



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cloward-Style Cervical Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cloward-Style Cervical Retractor

SKU RM-CSCR-17-1

\$75 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cloward-Style Cervical Retractor Frame	17 cm overall length, hinged arms, ratchet lock	\$75
Protective Sleeve / Box	Individual transit packaging	\$75

PRODUCT DESCRIPTION

The Cloward-Style Cervical Retractor is a foundational self-retaining surgical retractor frame engineered specifically to provide steady, hands-free soft-tissue displacement during anterior cervical spine procedures. Built with an articulated hinged architecture, this frame conforms dynamically to patient-specific neck contours, distributing retraction forces evenly while protecting delicate tracheoesophageal tissues.

Designed to eliminate manual assistant fatigue, this 17 cm Cloward-Style Cervical Retractor chassis features a hinged frame-only layout (\$75) engineered to accept interchangeable cervical blades. The finger-ring handle integrates a multi-position mechanical ratchet lock assembly that holds targeted tissue tension automatically and securely, preventing accidental slipping while minimizing hand and wrist fatigue during prolonged discectomy and fusion maneuvers.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the Cloward-Style Cervical Retractor non-sterile and fully autoclavable inside protective packaging, providing spine surgeons and neurosurgical teams with dependable durability, smooth ratchet action, and long-term clinical performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)