



RETRACTOR MAKER

PRODUCT BROCHURE OF

Humeral Head Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Humeral Head Retractor

SKU RM-HHR-ORTHO-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8 ¾" (22.23 cm) — Curved shank, angled ring blade	RM-HHR-ORTHO-1-1	\$84
9.0" (22.86 cm) — Curved shank, angled ring blade	RM-HHR-ORTHO-1-2	\$89
10.0" (25.40 cm) — Curved shank, angled ring blade	RM-HHR-ORTHO-1-3	\$94
11.0" (27.94 cm) — Curved shank, angled ring blade	RM-HHR-ORTHO-1-4	\$99

PRODUCT DESCRIPTION

The Humeral Head Retractor is a specialized orthopedic surgical instrument engineered specifically to provide reliable bone leverage and soft-tissue displacement during total and reverse shoulder arthroplasty. Designed with an anatomically curved shank and an angled ring-blade profile, this retractor cups the proximal humerus securely, levering it posteriorly away from the glenoid rim without slipping into surrounding neurovascular zones.

Built for complete procedural adaptability across variable patient sizes, the Humeral Head Retractor series is offered in four selectable overall tool lengths: 8 ¾ inches (\$84), 9 inches (\$89), 10 inches (\$94), and an extended 11 inches (\$99) for deep muscular or bariatric shoulder corridors. Its contoured blade edges distribute leverage counter-traction evenly across cortical bone, while the heavy-duty handle provides surgical assistants with sustained manual pulling power, reducing hand fatigue throughout prolonged glenoid reaming and component placement.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Humeral Head Retractor to withstand immense orthopaedic torque and repeated hospital autoclave cycles, providing shoulder reconstruction and sports medicine teams worldwide with unyielding clinical durability, corrosion resistance, and surgical precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)