



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hohmann Retractor Right Angled

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hohmann Retractor Right Angled

SKU RM-HRRA-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
9 ½" Overall Length — 24.13 cm length, 7" handle, 40 mm right-angle bend	RM-HRRA-1-1	\$84
9 ¾" Overall Length — 24.76 cm length, 7" handle, 18 mm right-angle bend	RM-HRRA-1-2	\$79

PRODUCT DESCRIPTION

The Hohmann Retractor Right Angled is a specialized orthopedic surgical instrument engineered specifically to provide direct, perpendicular tissue leverage and bone exposure during hip reconstruction and complex pelvic osteotomies. Designed with a true right-angled blade drop, this retractor allows surgical teams to retract deep muscle planes cleanly without crowding the operative corridor or exerting uneven pressure vectors.

Built to adapt to specific structural depths, the Hohmann Retractor Right Angled is offered in two distinct longitudinal size variants. The 9 ½-inch (24.13 cm) model (\$84) features a 7-inch handle configured to bend at a deep 40 mm right angle, ideal for profound pelvic and posterior hip exposures. The 9 ¾-inch (24.76 cm) model (\$79) features a 7-inch handle configured to bend at an 18 mm right angle for intermediate depths. Its contoured tip beak hooks securely over bone margins to prevent accidental slippage under heavy counter-traction.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Hohmann Retractor Right Angled to withstand heavy mechanical torque and repeated hospital autoclave cycles, providing orthopedic joint replacement and trauma teams worldwide with unyielding clinical durability and precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)