



RETRACTOR MAKER

PRODUCT BROCHURE OF

Casper Cervical Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron

Casper Cervical Retractor

SKU RM-CCR-115-1

\$199 USD



AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Casper Cervical Retractor Chassis	11.5 cm frame, 4-arm hinged transversal design, 85 mm spread	\$199
Protective Sleeve / Box	Individual transit packaging	\$199

PRODUCT DESCRIPTION

The Casper Cervical Retractor is an advanced self-retaining surgical retractor platform engineered specifically to provide reliable, multi-planar soft-tissue exposure during anterior cervical spine surgery. Designed with an articulated four-arm transversal chassis, this instrument conforms smoothly over the patient's neck, distributing retraction force evenly across wound margins while shielding vulnerable visceral structures.

Built to clinical framework standards, this Casper Cervical Retractor model features an overall frame length of 11.5 cm (\$199) capable of achieving a maximum opening spread of 85 mm. Its four hinged arms allow versatile intraoperative positioning, adapting to patient-specific lordosis while preventing localized muscle ischemia. The integrated mechanical ratchet-lock mechanism sustains continuous retraction tension automatically, eliminating the need for manual hand placement and providing uncompromised stability for cage insertion.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker manufactures the Casper Cervical Retractor to withstand repetitive hospital sterilization cycles, providing spine surgeons with exceptional mechanical rigidity, smooth ratchet control, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)