



RETRACTOR MAKER

PRODUCT BROCHURE OF

Tracheal Hook/Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Tracheal Hook/Retractor

SKU RM-THR-6IN-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard — 6 in, 15.2 cm	RM-THR-6IN-1-1	\$29
Large — 7 in, 18 cm	RM-THR-6IN-1-2	\$34
X-Large — 7.5 in, 19 cm	RM-THR-6IN-1-3	\$39

PRODUCT DESCRIPTION

The Tracheal Hook/Retractor is an essential emergency and elective airway surgical instrument engineered specifically to provide safe, non-slip elevation and stabilization of tracheal cartilages during tracheostomy procedures. Designed with an anatomically curved sharp hook, this instrument allows surgeons to engage and elevate the cricoid or upper tracheal rings cleanly, lifting the airway toward the operative incision while displacing overlying pretracheal fascia.

Built to withstand high-frequency autoclave sterilization cycles without degrading, this Tracheal Hook/Retractor configuration features an overall tool length of 6 inches (15.2 cm, \$29). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure lasting corrosion resistance, its working tip profile incorporates a single-ended design featuring a curved sharp single prong configuration that bites securely into dense perichondrium, preventing instrument displacement during critical tracheostomy incision and cannula insertion. The ergonomic flat handle delivers steady manual leverage for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Tracheal Hook/Retractor to uphold strict clinical standards, providing ENT suites, emergency departments, and critical care operating rooms worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)