



RETRACTOR MAKER

PRODUCT BROCHURE OF

Scoville Nerve Root Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Scoville Nerve Root Retractor

SKU RM-SNRR-235-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Scoville Nerve Root Retractor	23.5 cm overall length, angled contoured tip, slender shaft	\$29
Protective Sleeve / Box	Individual transit packaging	\$29

PRODUCT DESCRIPTION

The Scoville Nerve Root Retractor is an internationally acclaimed neurosurgical and spinal surgical retractor engineered specifically to provide safe, non-traumatic retraction and shielding of traversing spinal nerve roots and the thecal sac during posterior lumbar discectomy. Designed with an elongated slender body shaft, this retractor allows spine surgeons to work comfortably in deep interlaminar windows without obstructing microscopic visualization or crowding micro-instruments.

Classified as a specialized Premium Operating Room grade nerve root retractor, this dedicated Scoville Nerve Root Retractor configuration features an overall tool length of 23.5 cm (\$29). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to provide a lifetime warranty against rust and defects, its structural geometry incorporates a slender body shaft terminated with an angled, contoured working tip layout specifically engineered to cradle nerve fibers smoothly, preventing localized traction neuropraxia or dural lacerations. Its ergonomic flat handle delivers refined tactile feedback.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Scoville Nerve Root Retractor to withstand high-frequency autoclave sterilization, providing neurosurgical suites, spine centers, and orthopedic operating rooms worldwide with exceptional durability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)