



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kasden Double Ended delicate retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kasden Double Ended delicate retractor

SKU RM-KDEDR-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Kasden Double Ended delicate retractor	Dimensions / Details: 6 ½" (16.5 cm) overall length, dual delicate blades, aluminum handle	\$29
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging	\$29

PRODUCT DESCRIPTION

The Kasden Double Ended delicate retractor is a specialized handheld surgical instrument engineered specifically to provide atraumatic soft-tissue and wound margin retraction across delicate cosmetic, ophthalmic, and plastic reconstructive procedures. Designed with an ultra-fine double-ended architecture, this retractor allows surgical teams to transition smoothly between superficial cutaneous flap elevation and deeper soft-tissue retraction without changing tools.

Built for tight anatomical access and delicate wound boundary management, this Kasden Double Ended delicate retractor configuration features an overall tool length of 6 ½ inches (16.5 cm, \$29). Its blade geometry incorporates a double-ended delicate structure handcrafted with Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) that prevents localized bruising or tissue shearing. The handle configuration features a lightweight, anodized aluminum grip framework that provides tactile sensitivity, optimal balance, and reduced assistant hand strain during lengthy micro-dissections.

Manufactured to strict clinical tolerances with a reflection-free satin matte finish, Retractor Maker produces the Kasden Double Ended delicate retractor to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for plastic and ENT operating suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)