



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Cottle Retractor Hook 2 Prong

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



## Cottle Retractor Hook 2 Prong

**SKU** RM-CRH2P-1

**\$49 USD**

### AVAILABLE SIZES & PRICING

| SIZE / OPTION                               | SKU / DETAILS | PRICE (USD) |
|---------------------------------------------|---------------|-------------|
| Ball End Prongs — 6 ¾" length, 12 mm width  | RM-CRH2P-1-1  | \$49        |
| Sharp End Prongs — 6 ¾" length, 12 mm width | RM-CRH2P-1-2  | \$49        |

### PRODUCT DESCRIPTION

The Cottle Retractor Hook 2 Prong is a specialized handheld surgical instrument engineered specifically to provide precise, stable soft-tissue retraction during rhinoplasty, facial plastic, and reconstructive procedures. Designed with an agile dual-pronged rake architecture, this instrument allows surgeons to hook and mobilize skin flaps, columellar tissues, and cartilaginous margins with high tactile sensitivity.

Built for micro-anatomical clearance, this Cottle Retractor Hook 2 Prong model features an overall length of 6 ¾ inches (171 mm) paired with a 12 mm wide working framework (\$49). Surgeons can customize the tip geometry based on clinical demands: Ball End Prongs engineered to safely hold soft mucosa or delicate skin flaps without risking deep structural puncture, or Sharp End Prongs optimized to provide a reliable mechanical hold on dense facial tissues or cartilage segments. Its flat ergonomic handle ensures fatigue-free control.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Cottle Retractor Hook 2 Prong to withstand repetitive hospital sterilization cycles, providing plastic and ENT surgical suites with superior corrosion resistance, unyielding structural rigidity, and complete autoclave durability.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)