



RETRACTOR MAKER

PRODUCT BROCHURE OF

Tebbetts Spatula Breast Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Tebbetts Spatula Breast Retractor

SKU RM-TSBR-125-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
25 mm wide × 185 mm long — 12 ½" (31.5 cm) — Standard adult subpectoral pocket elevation and inframammary access	RM-TSBR-125-1-1	\$89
38 mm wide × 190 mm long — 12 ½" (31.5 cm) — Broad pectoral muscle displacement and large reconstructive pockets	RM-TSBR-125-1-2	\$99

PRODUCT DESCRIPTION

The Tebbetts Spatula Breast Retractor is an indispensable aesthetic and reconstructive plastic surgery retractor engineered specifically to provide broad, gentle tissue elevation during breast augmentation and mammaplasty. Designed with an elongated semi-malleable spatula blade, this instrument allows surgeons to adapt the blade curvature intraoperatively, following the natural convexities of the thoracic wall and rib cage without exerting localized ischemic pressure on breast parenchyma.

Built to withstand frequent autoclave sterilization cycles while providing excellent corrosion resistance, this Tebbetts Spatula Breast Retractor configuration features an overall tool length of 12 ½ inches (31.5 cm). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry layout incorporates a 1 × 6-inch single-ended semi-malleable spatula blade available across two separate width profiles: a 25 mm wide × 185 mm long blade footprint (\$89) for standard inframammary corridors, and an expanded 38 mm wide × 190 mm long blade footprint (\$99) for broad muscular elevation. The ergonomic flat handle ensures fatigue-free manual leverage for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Tebbetts Spatula Breast Retractor to withstand repetitive hospital sterilization cycles, providing cosmetic, plastic, and oncologic breast reconstructive suites worldwide with exceptional ductility, structural reliability, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)