



RETRACTOR MAKER

PRODUCT BROCHURE OF

Greater Trochanteric Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Greater Trochanteric Retractor

SKU RM-GTR-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Greater Trochanteric Retractor	14.875" (37.8 cm) length, 25 mm width × 1.5" depth blade — 1	\$49
Protective Sleeve / Box	Individual transit packaging — 1	\$49

PRODUCT DESCRIPTION

The Greater Trochanteric Retractor is a heavy-duty orthopedic surgical retractor engineered specifically to provide reliable bone leverage and soft-tissue displacement around the proximal femur and hip joint during total hip arthroplasty. Designed to hook over the lateral trochanteric ridge, this instrument allows surgical assistants to retract the gluteus medius, minimus, and vastus lateralis muscles away from the operative field without tearing tendon insertions.

Built for targeted exposure of the greater trochanter area, this Greater Trochanteric Retractor configuration features an extended overall tool length of 14.875 inches (37.8 cm, \$49) that keeps assistant hands well outside the surgical corridor. The working end incorporates a 25 mm wide blade profile with a 1.5-inch (3.8 cm) blade depth specifically contoured to seat securely over trochanteric bone margins without slipping. Its ergonomically designed handle maximizes grip control and minimizes user hand fatigue during prolonged joint reconstructions.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Greater Trochanteric Retractor to endure heavy orthopaedic torque and repeated hospital autoclave cycles, providing hip arthroplasty and trauma reconstruction teams worldwide with unyielding clinical durability and performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)