



RETRACTOR MAKER

PRODUCT BROCHURE OF

Gorney Rake Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Gorney Rake Retractor

SKU RM-GRR-3IN-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
3 in (7.5 cm), 38 mm wide — Standard	RM-GRR-3IN-1-1	\$49
4 1/4 in (10.8 cm), 38 mm wide	RM-GRR-3IN-1-2	\$59

PRODUCT DESCRIPTION

The Gorney Rake Retractor is an internationally acclaimed aesthetic and plastic surgery retractor engineered specifically to provide firm, non-slip skin flap elevation during rhytidectomy, neck lift, and reconstructive soft-tissue procedures. Designed with an expansive 4-prong rake head, this instrument allows surgeons to apply uniform manual counter-traction across large tissue planes, preventing focal tearing while maximizing operative visualization of underlying fascial layers.

Built for larger tissues requiring careful retraction, this Gorney Rake Retractor configuration features an overall blade length of 3 inches (7.5 cm, \$49) paired with an expansive 38 mm blade width that offers an extended surface contact area. Its working end incorporates a 4-prong configuration layout with finely sharpened tips engineered to grip dense dermal and subcutaneous tissue securely without slipping. The contoured handle provides surgical assistants with fatigue-free manual leverage throughout prolonged aesthetic dissections.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Gorney Rake Retractor to withstand repetitive hospital sterilization cycles, providing cosmetic, plastic, and reconstructive surgery centers worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)