



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cervical Tissue Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Cervical Tissue Retractor

SKU RM-CTR-13-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cervical Tissue Retractor	13 cm overall length, 1×1 blunt prong tip	\$79
Protective Sleeve / Box	Individual transit packaging	\$79

PRODUCT DESCRIPTION

The Cervical Tissue Retractor is a specialized handheld surgical instrument engineered specifically to provide localized, non-traumatic soft-tissue displacement during delicate anterior cervical spine exposures. Designed to maneuver within tight anatomical workspaces, this retractor allows surgical teams to gently mobilize fascia, muscular margins, and fatty tissue away from the longus colli muscles.

Built to optimize control in localized operative fields, this Cervical Tissue Retractor model features an overall length of 13 cm (\$79) paired with a 1×1 blunt prong configuration. The single blunt tip is carefully structured to hook and stabilize tissue margins effectively, preventing unnecessary patient tissue trauma or accidental puncture of the esophagus and carotid sheath. Its ergonomic flat handle provides tactile sensitivity and non-slip control during prolonged cervical procedures.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cervical Tissue Retractor to withstand repetitive hospital sterilization cycles, providing spinal and neurosurgical operating theaters with exceptional durability, corrosion resistance, and surgical reliability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)