



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sistrunk Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sistrunk Retractor

SKU RM-STR-153-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Sistrunk Retractor	15.3 cm overall length, 3/8" (9.5 mm) double-ended blades, flat handle	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Sistrunk Retractor is an indispensable head, neck, and endocrine surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement during thyroglossal duct cyst excision (Sistrunk procedure) and cervical dissections. Designed with a double-ended blade architecture, this instrument allows surgeons to elevate skin edges, cervical fascia, and strap muscles cleanly, establishing an unhindered visual corridor to the hyoid bone without compressing delicate anterior neck structures.

Built to deliver maximum durability and corrosion resistance backed by an absolute lifetime warranty against rust, this Sistrunk Retractor configuration features an overall tool length of 15.3 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade width incorporates a 3/8-inch (approx. 9.5 mm) wide blade layout contoured with smooth, rounded margins to protect delicate tissue flaps from tearing. The structural geometry integrates a double-ended blade chassis into a flat central grip handle that provides operating assistants with fatigue-free manual leverage.

Supplied as a non-sterile, re-sterilizable operating room device with an anti-glare satin matte finish, Retractor Maker manufactures the Sistrunk Retractor to withstand frequent sterilization runs, providing ENT operating suites, endocrine centers, and plastic surgery clinics worldwide with unyielding structural reliability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)