



RETRACTOR MAKER

PRODUCT BROCHURE OF

Harrington Splanchnic Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Harrington Splanchnic Retractor

SKU RM-HSR-SPL-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
127 mm long × 38 mm wide — 30.5 cm or 31.8 cm	RM-HSR-SPL-1-1	\$69
127 mm long × 63 mm wide — 30.5 cm or 31.8 cm	RM-HSR-SPL-1-2	\$69
190 mm long × 63 mm wide — 30.5 cm or 31.8 cm	RM-HSR-SPL-1-3	\$74

PRODUCT DESCRIPTION

The Harrington Splanchnic Retractor is a premier visceral surgical retractor engineered specifically to provide broad, atraumatic displacement of parenchymal abdominal organs and deep pelvic structures during exploratory laparotomy and oncologic surgery. Designed with an iconic sweetheart heart-shaped distal paddle, this instrument distributes manual counter-traction evenly across a broad surface area, protecting fragile visceral structures such as the liver, spleen, and intestinal loops from localized pressure trauma.

Built for variable operative depths, this Harrington Splanchnic Retractor series is available in two overall lengths: 30.5 cm (\$69) and 31.8 cm (\$74). Surgeons can customize the blade footprint across three distinct dimensions: 127 mm long × 38 mm wide for narrow fields, 127 mm long × 63 mm wide for standard laparotomies, and an extra-deep 190 mm long × 63 mm wide profile for bariatric or retroperitoneal procedures. Its ergonomic handle provides surgical assistants with steady, non-slip manual leverage.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Harrington Splanchnic Retractor to withstand heavy mechanical loads and repeated hospital autoclave cycles, delivering exceptional corrosion resistance, structural rigidity, and complete clinical durability for general and gynecological surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)