



RETRACTOR MAKER

PRODUCT BROCHURE OF

Deep Rake Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Deep Rake Retractor

SKU RM-DRR-1

\$68 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
38 mm Depth — 57 mm wide × 38 mm deep — Intermediate subcutaneous and muscular retraction	RM-DRR-1-1	\$68
51 mm Depth — 57 mm wide × 51 mm deep — Standard abdominal wall and orthopedic joint exposure	RM-DRR-1-2	\$72
63 mm Depth — 57 mm wide × 63 mm deep — Deep pelvic floor and paraspinal muscle displacement	RM-DRR-1-3	\$76
76 mm Depth — 57 mm wide × 76 mm deep — Maximum-depth bariatric and heavy muscular leverage	RM-DRR-1-4	\$80

PRODUCT DESCRIPTION

The Deep Rake Retractor is an indispensable surgical instrument engineered specifically to provide powerful, controlled tissue displacement during deep abdominal, orthopedic, and pelvic operations. Designed with a robust 4-prong blunt rake architecture, this retractor allows surgical teams to mobilize thick subcutaneous fat layers, heavy muscle masses, and fibrous fascial planes without tearing tissue.

Built for versatile deep-cavity clearance, this Deep Rake Retractor series features an extended overall length of 27 cm that keeps the assistant's hands safely outside the surgical corridor. The instrument incorporates a 57 mm wide rake head available across four selectable blade depths: 57 mm × 38 mm (\$68), 57 mm × 51 mm (\$72), 57 mm × 63 mm (\$76), and 57 mm × 76 mm (\$80). Its rounded blunt prongs distribute traction forces

evenly, while the ergonomic teardrop ring handle provides secure, fatigue-free manual leverage.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Deep Rake Retractor to strict clinical standards, delivering high tensile strength, exceptional corrosion resistance, and complete autoclave durability for general, orthopedic, and trauma surgical suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)