



RETRACTOR MAKER

PRODUCT BROCHURE OF

Parsonnet Epicardial Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Parsonnet Epicardial Retractor

SKU RM-PER-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Extra Small	RM-PER-1-1	\$44
Small	RM-PER-1-2	\$49
Medium	RM-PER-1-3	\$54
Large	RM-PER-1-4	\$59

PRODUCT DESCRIPTION

The Parsonnet Epicardial Retractor is a specialized cardiovascular surgical retractor engineered specifically to provide gentle, hands-free displacement of epicardial fat and myocardial tissue overlying coronary vessels. Designed with a flexible, spring-loaded wire framework, this instrument isolates target coronary arteries cleanly, establishing a stable, bloodless operative corridor during delicate arteriotomies and microvascular anastomoses.

Built to accommodate varied patient anatomical boundaries, this Parsonnet Epicardial Retractor line keeps a uniform spring design while offering four selectable sizes: Extra Small at 2.5 cm (\$44) for fine distal vessels, Small at 3.0 cm (\$49), Medium at 3.8 cm (\$54) for standard coronary exposures, and Large at 4.8 cm (\$59) for deep intramyocardial vessels. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its contoured, polished wire prongs embed securely into the epicardium without puncturing myocardium or compressing coronary lumina.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Parsonnet Epicardial Retractor to withstand repetitive hospital sterilization cycles, delivering high elastic memory, corrosion resistance, and complete autoclave durability for cardiac surgical suites and congenital heart institutes worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)