



RETRACTOR MAKER

PRODUCT BROCHURE OF

Krasky Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Krasky Retractor

SKU RM-KRK-28-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Krasky Retractor	Dimensions / Details: 28 cm overall length, 50 × 25 mm blade, 6 sharp prongs · Quantity Included: 1	\$99
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$99

PRODUCT DESCRIPTION

The Krasky Retractor is an indispensable deep-cavity surgical instrument engineered specifically to provide dependable, non-slip muscular displacement and wound margin exposure during major abdominal and orthopedic procedures. Designed with an elongated curved blade, this retractor allows surgical assistants to apply steady counter-traction across heavy tissue masses without slipping into the operative field.

Built for heavy surgical counter-traction, this dedicated Krasky Retractor configuration features an overall tool length of 28 cm (\$99). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its working end incorporates a 50 mm long by 25 mm wide curved blade geometry equipped with 6 sharp rake-style prongs that bite securely into dense fibrous tissues and thick muscle planes. The ergonomic handle provides operating assistants with superior manual leverage, reducing hand fatigue during multi-hour surgical exposures.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Krasky Retractor to withstand repetitive hospital autoclave cycles, providing general, orthopedic, and trauma surgical teams worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)