



RETRACTOR MAKER

PRODUCT BROCHURE OF

Henly Retractor Center Blade

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Henly Retractor Center Blade

SKU RM-HRCB-1

\$48 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
19 mm long × 15 mm wide — Standard Henly stem mount	RM-HRCB-1-1	\$48
25 mm long × 15 mm wide — Standard Henly stem mount	RM-HRCB-1-2	\$50
32 mm long × 15 mm wide — Standard Henly stem mount	RM-HRCB-1-3	\$52
76 mm long × 15 mm wide — Standard Henly stem mount	RM-HRCB-1-4	\$54

PRODUCT DESCRIPTION

The Henly Retractor Center Blade is a specialized modular surgical component engineered specifically to provide secure, non-traumatic center-wall retraction when mounted onto a compatible Henly self-retaining retractor frame. Designed for deep vascular and soft-tissue exposure, this blade transforms standard two-way lateral spreading into a stable three-way operative corridor, shielding critical structures from accidental injury.

Built for complete procedural adaptability, this Henly Retractor Center Blade features a standardized mounting stem that locks firmly into the retractor assembly crossbar. Surgeons can select from four distinct blade dimensions: 19 mm long × 15 mm wide (\$48), 25 mm long × 15 mm wide (\$50), 32 mm long × 15 mm wide (\$52), and an extended 76 mm long × 15 mm wide (\$54) profile for deep femoral and retroperitoneal corridors. Its curved lip cradles tissue margins securely without causing ischemic pressure points.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Henly Retractor Center Blade to withstand repetitive hospital autoclave

cycles, providing vascular, general, and orthopedic surgical suites worldwide with unyielding structural rigidity, high corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)