



RETRACTOR MAKER

PRODUCT BROCHURE OF

Millard Cottle Retractor with Finger Ring

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Millard Cottle Retractor with Finger Ring

SKU RM-MCRFR-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Millard Cottle Retractor with Finger Ring	5.5 cm length, 5 × 12 mm sharp 2-prong blade, finger ring	\$24
Protective Sleeve / Box	Individual transit packaging	\$24

PRODUCT DESCRIPTION

The Millard Cottle Retractor with Finger Ring is a specialized rhinoplasty surgical instrument engineered specifically to provide superior manual dexterity and stable alar retraction within confined intranasal spaces. Designed with an ultra-compact chassis, this retractor allows the operating surgeon to maintain continuous counter-traction on the nostril rim and lower lateral cartilage while keeping the operative corridor unobstructed for delicate trimming and suturing.

Built for high-maneuverability facial micro-exposure, this Millard Cottle Retractor with Finger Ring model features an overall length of 5.5 cm (\$24). Its distal working end incorporates a 5 mm × 12 mm sharp double-prong blade that engages dense fibrous tissues and vestibular skin securely without slipping. The integrated finger ring at the base of the handle stem provides exceptional tactile feedback, stabilizes the assistant or surgeon's hand, and prevents accidental slippage under wet operative conditions.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Millard Cottle Retractor with Finger Ring to withstand repetitive hospital sterilization cycles, providing plastic, aesthetic, and ENT surgeons with unyielding clinical durability and micro-dissection control.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)