



RETRACTOR MAKER

PRODUCT BROCHURE OF

Martin Arm Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Martin Arm Retractor

SKU RM-MAR-1

\$799 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Coupler	Clamps securely onto the OR table side rail while holding and locking the L Shape Bar in place
1	L Shape Bar	Inserts into the Coupler to serve as the primary vertical and horizontal structural support framework
1	Articulating Arm	Connects to the L Shape Bar to provide multi-axis movement and locked positioning for retractor instruments over the surgical site

PRODUCT DESCRIPTION

The Martin Arm Retractor is an indispensable table-mounted mechanical holding device designed to eliminate assistant strain and guarantee tremor-free instrument stability during open and laparoscopic interventions. By anchoring directly onto standard operating table accessory side rails, this rigid platform serves as a dependable static holding arm for heavy tissue retractors, endoscopes, and optical camera heads.

Built for seamless multi-planar positioning, the Martin Arm Retractor incorporates a universal table rail coupler, a heavy-duty vertical and horizontal L-shaped support bar, and an advanced articulating friction arm. With 360-degree rotational flexibility across multiple locking joints, the surgeon can position retraction blades or diagnostic scopes at any precise angle and lock them securely into place using an ergonomic central tightening handle.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a glare-reducing satin finish, Retractor Maker supplies the Martin Arm Retractor non-sterile and fully autoclavable inside a protective storage box, delivering rugged clinical longevity to surgical facilities worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)