



RETRACTOR MAKER

PRODUCT BROCHURE OF

Shoulder Arthroplasty Retractor Glenoid Exposure

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Shoulder Arthroplasty Retractor Glenoid Exposure

SKU RM-SAR-GE-775-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Shoulder Arthroplasty Retractor Glenoid Exposure	7 ¾" (19.68 cm) length, 40 mm wide blade, ergonomic T-handle	\$79
Protective Sleeve / Box	Individual transit packaging	\$79

PRODUCT DESCRIPTION

The Shoulder Arthroplasty Retractor Glenoid Exposure is a specialized orthopedic surgical instrument engineered specifically to provide powerful bone leverage and circumferential soft-tissue displacement during total and reverse shoulder reconstructions. Designed to seat behind the posterior glenoid rim and humeral neck, this retractor allows surgical teams to displace the proximal humerus posteriorly, establishing an expansive visual corridor directly into the glenoid vault for reaming and component seating.

Built to withstand high-frequency autoclave sterilization protocols, this Shoulder Arthroplasty Retractor Glenoid Exposure configuration features an overall tool length of 7 ¾ inches (19.68 cm, \$79). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade width incorporates an expansive 40 mm wide blade layout contoured to cradle heavy capsular and muscular tissues without crushing bone cortices. Its integrated ergonomic T-handle is engineered to promote a powerful grip and enhance mechanical leverage for the user, significantly reducing assistant forearm fatigue during prolonged glenoid preparation.

Provided as a non-sterile, re-sterilizable operating room device with an anti-glare satin matte finish, Retractor Maker manufactures the Shoulder Arthroplasty Retractor Glenoid Exposure to withstand heavy orthopedic torque, delivering exceptional corrosion resistance, structural rigidity, and complete autoclave durability for shoulder reconstruction suites and sports medicine centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)