



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hohmann Retractors Square Blade

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hohmann Retractors Square Blade

SKU RM-HRSB-1

\$74 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8 ½" (21.59 cm) — Square blade profile, extra-long handle	RM-HRSB-1-1	\$74
9 ½" (24.13 cm) — Square blade profile, extra-long handle	RM-HRSB-1-2	\$79

PRODUCT DESCRIPTION

Hohmann Retractors Square Blade are specialized orthopedic surgical instruments engineered specifically to provide broad-surface bone contact and powerful muscle displacement during joint reconstruction and fracture fixation. Designed with a wide, planar square blade face, these retractors seat securely over long bone cortices without causing localized pressure fractures or gouging fragile osteoporotic bone.

Built for superior mechanical advantage, the Hohmann Retractors Square Blade series features an extra-long handle layout across all versions, allowing surgical assistants to apply sustained counter-traction while keeping hands outside the line of sight. Surgeons can choose between two distinct overall lengths: 8 ½ inches (21.59 cm, \$74) and 9 ½ inches (24.13 cm, \$79). The square blade profile hooks firmly behind cortical margins, preventing tissue creep and slippage during vigorous bone reaming and plating.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces Hohmann Retractors Square Blade instruments to withstand heavy orthopaedic torque and repeated hospital autoclave cycles, providing trauma and joint replacement surgical teams worldwide with exceptional clinical durability and performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)