



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ranawat Knee Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Ranawat Knee Retractor

SKU RM-RKR-95-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Ranawat Knee Retractor	9 ½" (24.13 cm) length, 17 × 34 mm 90° blade (Blunt or Sharp)	\$79
Protective Sleeve / Box	Individual transit packaging	\$79

PRODUCT DESCRIPTION

The Ranawat Knee Retractor is a specialized orthopedic surgical instrument engineered specifically to provide dependable bone leverage and soft-tissue displacement during total knee arthroplasty (TKA) and complex knee reconstructions. Designed to drop behind the tibial plateau and collateral ligaments, this retractor allows surgical teams to hold back heavy capsular and ligamentous structures cleanly away from oscillating bone saws and reamers.

Built with a highly durable and corrosion-resistant chassis, this Ranawat Knee Retractor configuration features an overall tool length of 9 ½ inches (24.13 cm, \$79). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its working tip dimensions measure 17 mm wide × 34 mm long. The tip profile geometry incorporates an integrated 90-degree angled working end available in two distinct typologies selectable from the option menu: a Blunt tip engineered for gentle, atraumatic retraction of vascularly sensitive soft tissues, or a Sharp tip engineered to provide a highly secure, non-slip purchase on tough ligamentous structures and cortical bone. The heavy-duty handle delivers fatigue-free manual leverage for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Ranawat Knee Retractor to withstand heavy orthopedic torque and repeated hospital autoclave cycles, providing knee arthroplasty centers and orthopedic trauma surgical suites worldwide with exceptional mechanical stability, corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)