



RETRACTOR MAKER

PRODUCT BROCHURE OF

Acetabular Retractor Narrow Double Prong

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Acetabular Retractor Narrow Double Prong

SKU RM-ARNDP-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Narrow — 32 mm blade	RM-ARNDP-1-1	\$99
Standard — 44 mm blade	RM-ARNDP-1-2	\$99

PRODUCT DESCRIPTION

The Acetabular Retractor Narrow Double Prong is a heavy-duty orthopedic surgical retractor engineered to provide exceptional visual access to the acetabulum during total hip replacement and hip reconstruction surgery. By anchoring its dual prongs securely against the acetabular margin or ischium, this instrument allows surgical teams to retract the femur and surrounding musculature with powerful mechanical leverage.

Designed for robust intraoperative stability, this Acetabular Retractor Narrow Double Prong configuration features an extended 12 ½-inch (31.75 cm) shaft that keeps the assistant's hands outside the direct visual field. The instrument is available in two specialized blade widths—32 mm and 44 mm—featuring a curved neck and narrow dual prongs that engage bone firmly without slipping or traumatizing adjacent neurovascular bundles.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a high-tensile core and anti-glare satin finish, Retractor Maker manufactures the Acetabular Retractor Narrow Double Prong to endure heavy orthopaedic torque and repeated steam sterilization cycles, providing joint reconstruction teams with long-lasting clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)