



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ferreira Breast Augmentation Retractor - Fiber Optics

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Ferreira Breast Augmentation Retractor - Fiber Optics

SKU RM-FBAR-FO-1

\$129 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Ferreira Breast Augmentation Retractor - Fiber Optics	5 1/8" length, 25 x 171 mm curved blade, fiber-optic bundle
1	Universal Light Cable Fitting	Connects to standard OR light cables
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Ferreira Breast Augmentation Retractor - Fiber Optics is a premier illuminated aesthetic surgery instrument engineered specifically to provide superior pocket visualization and gentle soft-tissue elevation during mammoplasty. By channeling high-intensity fiber-optic illumination along the undersurface of its blade, this instrument eliminates dark recesses in deep subpectoral pockets, allowing surgeons to verify hemostasis and position implants with absolute confidence.

Built for maximum cavity clearance and safety, this specialized Ferreira Breast Augmentation Retractor - Fiber Optics model features an overall tool length of 5 1/8 inches (13.0 cm, \$129) paired with a 25.0 mm wide by 171.0 mm long blade. Its curved and tapered layout ensures gentle skin deflection and soft tissue manipulation, preventing shear stresses on cutaneous borders. The integrated fiber-optic light bundle features a universal fitting compatible with standard operating room light source cables (Storz, Olympus, Wolf, ACMI).

Crafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-transmission optical glass fibers with an anti-glare satin finish, Retractor Maker supplies the Ferreira Breast Augmentation Retractor - Fiber Optics non-sterile and fully autoclavable inside protective packaging, providing cosmetic, plastic, and oncologic breast surgeons with high durability and optical clarity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)