



RETRACTOR MAKER

PRODUCT BROCHURE OF

Czerny Retractor Double-Ended

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Czerny Retractor Double-Ended

SKU RM-CRDE-CZ-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Czerny Retractor Double-Ended	17.5 cm overall length, 38 × 22 mm dual-ended blades	\$44
Protective Sleeve / Box	Individual transit packaging	\$44

PRODUCT DESCRIPTION

The Czerny Retractor Double-Ended is an adaptable handheld surgical retractor engineered specifically to provide dependable, multi-depth tissue displacement across abdominal, pelvic, and general surgical procedures. Built with an S-shaped dual-ended architecture, this instrument allows surgical teams to transition smoothly between superficial wound margin retraction and deeper fascial holding without having to switch instruments during active dissection.

Designed to maximize clinical versatility, this Czerny Retractor Double-Ended configuration features an overall length of 17.5 cm (\$44) paired with 38 mm × 22 mm blade jaws. Its dual-ended design incorporates two distinct blade geometries to seamlessly handle differing tissue thicknesses and localized depth requirements. The gently cupped, lip-ended blade on one end hooks securely behind muscular layers, while the opposing right-angled flat blade maintains clean visual margins over cutaneous and subcutaneous planes.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Czerny Retractor Double-Ended to withstand repetitive hospital autoclave cycles, providing operating theaters worldwide with high corrosion resistance, unyielding structural rigidity, and dependable clinical performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)