



RETRACTOR MAKER

PRODUCT BROCHURE OF

Finochietto-Infant Rib Spreader

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Finochietto-Infant Rib Spreader

SKU RM-FIRS-SERIES-1

\$319 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Delicate Horizontal Track Frame	Compact slide bar layout optimized for infantile thoracotomy
2 (Fixed Pair)	Fixed Angled Blades	14×15 mm or 15×15 mm fenestrated lateral modules
1	Turning Screw Lock	Manual safety locking bolt that fixes retraction and limits fatigue
2	Ergonomic Tapered Forceps Arms	Lightweight structured handles built with completely smooth outer edges
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Finochietto-Infant Rib Spreader is an essential neonatal and infantile thoracic surgical retractor engineered specifically to provide controlled, micro-calibrated rib distraction during congenital cardiothoracic procedures. Designed with a delicate horizontal slide track, this instrument allows pediatric surgeons to separate tiny infant ribs with exceptional gentleness, avoiding costochondral junction tearing while maintaining an open visual corridor into the infant mediastinum.

Built to accommodate delicate infantile anatomical corridors, this Finochietto-Infant Rib Spreader series offers up to 75 mm of structural widening capacity across two specialized blade configurations: Blade Option A (\$319) features 14 mm × 15 mm fixed fenestrated blades, and Blade Option B (\$349) incorporates 15 mm × 15 mm fixed fenestrated blades. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the compact lateral modules cradle infant ribs securely without slipping. Its manual turning screw lock safety bolt fixes retraction securely to eliminate user fatigue, while ergonomic tapered forceps arms feature smooth outer

edges to protect adjacent soft tissues.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Finochietto-Infant Rib Spreader to withstand repetitive hospital sterilization cycles, providing neonatal surgical suites, congenital heart centers, and pediatric trauma theaters worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)