



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ring Breast Retractor With Fiber Optic Illumination

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Ring Breast Retractor With Fiber Optic Illumination

SKU RM-RBR-FO-339-1

\$339 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Ring Breast Retractor With Fiber Optic Illumination	Handheld breast retractor, ring blade, integrated light bundle
1	Universal Light Cable Adapter	Connects to standard OR fiber-optic cables
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Ring Breast Retractor With Fiber Optic Illumination is a premier illuminated aesthetic and reconstructive plastic surgery retractor engineered specifically to combine broad tissue elevation with direct, shadowless optical illumination inside deep breast pockets. By channeling high-intensity fiber-optic light directly along the blade undersurface, this instrument eliminates optical shadows in submuscular and subglandular spaces, allowing surgeons to verify complete hemostasis and precise pocket borders.

Built for superior surgical access during breast augmentation and reconstruction, this Ring Breast Retractor With Fiber Optic Illumination model (\$339) features an anatomically curved ring-shaped blade face that cradles breast parenchyma and the pectoralis major muscle atraumatically. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) paired with high-transmission optical glass fibers, the instrument incorporates an integrated universal connector compatible with standard OR light cables (Storz, Olympus, Wolf, ACMI). Its ergonomic handle ensures fatigue-free manual leverage for operating assistants.

Forged with an anti-glare satin matte finish on all exposed metallic surfaces, Retractor Maker manufactures the Ring Breast Retractor With Fiber Optic Illumination to withstand repetitive hospital autoclave cycles, providing plastic, aesthetic, and oncologic breast surgery centers worldwide with exceptional optical transmission, unyielding structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)