



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hoyos Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hoyos Retractor System

SKU RM-HRS-1

\$1,999 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Hoyos Distal Retractor Short	Short-length distal retraction blade option for access closer to incisions
1	Hoyos Distal Retractor Long	Long-length distal retraction blade option for deeper tissue reach
1	Hoyos Proximal Retractor Long	Long-length proximal retraction blade component for proximal boundary elevation
1	Remora Retractor	Specialized traction retractor blade designed to reduce slippage
1	Interlocking Support Rod	Connects retractor arms to the main overhead framework
1	Table Clamp for Hoyos Retractor System	Anchors the vertical post to the operating table rail
1	Hoyos Horizontal Extension Bar	Extends the framework horizontally over the surgical field
1	Hoyos Vertical Post	Primary vertical column supporting the mounting frame

PRODUCT DESCRIPTION

The Hoyos Retractor System is a premier table-mounted surgical retractor platform engineered to meet the demanding soft-tissue exposure requirements of advanced aesthetic and reconstructive plastic surgery. By

mounting rigidly to the operating table side rail, this system eliminates manual retraction fatigue, allowing surgical teams to maintain a continuous, symmetrical tissue flap elevation during extensive abdominal and body-contouring operations.

Designed specifically for modern aesthetic procedures, the Hoyos Retractor System features an 8-piece modular framework comprising a table clamp, vertical support post, horizontal extension bar, and an interlocking support rod. The retractor assortment includes a short distal retractor, a long distal retractor, a long proximal retractor, and the specialized anti-slip Remora retractor blade engineered to maintain secure grip on subcutaneous adipose tissue and fascial planes without tissue necrosis.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-strength metal alloys with an anti-glare satin finish, Retractor Maker supplies the complete Hoyos Retractor System non-sterile and fully autoclavable inside a heavy-duty aluminum sterilization container, ensuring rapid turnover and maximum procedural precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)