



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

FK Retractor (Feyh-Kastenbauer)

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# FK Retractor (Feyh-Kastenbauer)

**SKU** RM-FK-1

**\$4,999 USD**

## COMPONENT LIST

| QUANTITY | COMPONENT                                          | DIMENSIONS / DESCRIPTION              |
|----------|----------------------------------------------------|---------------------------------------|
| 1        | FK-WO TORS Basic Frame with Adapter                | Main structural frame assembly        |
| 1        | Adapter for FK-WO TORS Basic Frame                 | Mounting connection adapter           |
| 1        | Adapter for Tumor Grasping Forceps or Suction Tube | Auxiliary instrument mounting adapter |
| 1        | Light Clip                                         | Lateral cable connection              |
| 2        | Cheek Retractors                                   | Curved profile                        |
| 1        | Suction Tube for Smoke                             | 4.0 × 170 mm                          |
| 10 pcs   | Teeth Protector                                    | Latex-free, non-autoclavable          |
| 1        | Mandible Blade                                     | 11 cm                                 |
| 1        | Laryngeal Blade                                    | 17 cm, concave                        |
| 1        | Tongue Blade                                       | Curved                                |
| 1        | Tongue Blade                                       | Curved, right side open               |
| 1        | Tongue Blade                                       | Curved, left side open                |
| 1        | TORS Blade according to Weinstein-O'Malley         | Standard TORS profile                 |

| QUANTITY | COMPONENT               | DIMENSIONS / DESCRIPTION          |
|----------|-------------------------|-----------------------------------|
| 2        | Small TORS Blade, Left  | Compact left-sided configuration  |
| 2        | Small TORS Blade, Right | Compact right-sided configuration |
| 1        | Sterilization Box       | Formed tray for FK-WO TORS set    |

## PRODUCT DESCRIPTION

The FK Retractor (Feyh-Kastenbauer) is an internationally recognized trans-oral operative platform engineered to provide optimal visualization and surgical access to the tongue base, hypopharynx, supraglottis, and larynx. Optimized for Trans-Oral Robotic Surgery (TORS) as well as endolaryngeal laser microsurgery, this system establishes an expansive, unyielding visual corridor that accommodates robotic arms and microscopic visualization without tissue obstruction.

Engineered around the FK-WO TORS basic frame, this complete FK Retractor (Feyh-Kastenbauer) set features specialized mounting adapters for grasping forceps or suction tubes, light cable clips, curved cheek retractors, a dedicated 4.0 × 170 mm smoke suction tube, and 10 latex-free dental protectors. The blade assembly provides unmatched clinical flexibility, incorporating an 11 cm mandible blade, a 17 cm concave laryngeal blade, standard and open-sided curved tongue blades, and bilateral Weinstein-O'Malley TORS blades.

Supplied in a custom-formed heavy-duty sterilization box, Retractor Maker manufactures the FK Retractor (Feyh-Kastenbauer) from high-grade biocompatible stainless steel and specialized alloys, delivering dependable stability and surgeon-focused ergonomics directly to otolaryngology departments worldwide.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)