



RETRACTOR MAKER

PRODUCT BROCHURE OF

Urology Lobster Expanded Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Urology Lobster Expanded Retractor

SKU RM-ULER-1

\$3,499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
8	Lobster Retractor Clamp	—
1	Lobster Table Clamp with Vertical Post & Frame Coupler	—
1	Lobster Hinged Extension Arm	—
1	Lobster Frame	—
2	Lobster Mayo Swivel Retractor	2" × 2"
2	Lobster Mayo Swivel Retractor	2¾" × 2"
1	Lobster Deaver Swivel Retractor	2¾" × 4½"
1	Lobster Deaver Swivel Retractor	1½" × 6"
1	Lobster Harrington Sweetheart Retractor	2" × 5"
1	Lobster Harrington Retractor	2.5" × 6"
1	Lobster Malleable Swivel Retractor (Perforated)	2" × 6½"
1	Lobster Malleable Swivel Retractor (Perforated)	3" × 6½"

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Lobster Marshall Malleable Bladder Retractor, Swivel	3" × 6"
1	Lobster Fence Swivel Retractor	4" × 6"
1	Lobster Splanchnic Swivel Retractor	2" × 4½"
1	Sterilization Box	—

PRODUCT DESCRIPTION

The Urology Lobster Expanded Retractor is a specialized table-mounted self-retaining retractor system designed to overcome the visibility and mechanical challenges of deep urological surgery. By mounting securely to standard operating table rails with an articulated post and coupler, this system establishes a stable perimeter that frees surgical personnel from physical retraction while delivering stationary visualization into deep pelvic cavities.

Engineered specifically for extended pelvic and urological fields, the Urology Lobster Expanded Retractor configuration integrates a structural Lobster frame, an articulated hinged extension arm for reach, and eight multi-position retractor clamps. The 11-piece blade set features specialized urology retraction tooling, including a 3" × 6" Marshall malleable bladder swivel retractor, perforated malleable strips, Harrington Sweetheart blades, Deaver, Mayo, Fence, and visceral Splanchnic retractors that adjust to deep anatomical recesses.

Built to clinical standards from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-strength alloys with an anti-glare satin finish, Retractor Maker supplies the complete Urology Lobster Expanded Retractor in an organized, heavy-duty sterilization box validated for repeated steam autoclaving.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)