



RETRACTOR MAKER

PRODUCT BROCHURE OF

Greenberg Simplicity Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Greenberg Simplicity Retractor

SKU RM-GSR-1

\$2,499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Greenberg Primary Bar	6" (15 cm) — single-point rigid anchor to skull clamp post
1	Greenberg Floating Secondary Bar	12" (30 cm) — focused secondary track for mounting flex arms
4	Greenberg Long Retractor Arms	9.5" (24 cm) — multi-axis friction flex arms for precise blade hold
1	Greenberg Hand Rest	10" (25 cm) — adjustable forearm support platform
1	Greenberg Pattie Tray	3" x 6¾" (7.5 x 16.8 cm) — mounts on frame to store neuro-patties and tools
1	Compact Sterilization Box	Custom-molded layout — compact, perforated organizer and sterilizing box
1	Malleable Tapered Retractor Blade	9" x ¼" x ¾" — extra-wide blade for broad cortical protection
1	Malleable Tapered Retractor Blade	9" x 5/32" x ⅝" — medium-wide blade for lobar retraction
1	Malleable Tapered Retractor Blade	9" x ⅛" x ½" — standard blade for deep interhemispheric access
1	Malleable Tapered Retractor Blade	9" x 3/32" x ⅜" — fine-tapered blade for narrow cisternal pathways

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Malleable Tapered Retractor Blade	9" × 1/16" × 1/4" — micro-blade for delicate nerve/vessel separation

PRODUCT DESCRIPTION

The Greenberg Simplicity Retractor is a streamlined, space-efficient neurosurgical retraction system designed for rapid setup and unobstructed access during targeted cranial procedures. Ideal for procedures requiring minimal operative hardware clutter around the surgical field, this system delivers uncompromising holding power while maintaining a lightweight, compact footprint.

Optimized for small craniotomies and focused exposures, the Greenberg Simplicity Retractor combines one 6-inch primary anchor bar, one 12-inch floating secondary bar, four 9.5-inch multi-axis friction flex arms, and a 10-inch ergonomic handrest. The set includes a frame-mountable pattie tray alongside five double-ended malleable tapered retractor blades (from 1.6 mm micro-tips to 19 mm broad blades) to address precise cortical and deep cisternal holding requirements.

Manufactured entirely from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish, Retractor Maker supplies the Greenberg Simplicity Retractor non-sterile and fully autoclavable inside a custom-molded, compact perforated sterilization box built for repeated OR cycles.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)