



RETRACTOR MAKER

PRODUCT BROCHURE OF

Anderson Bear Claw Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Anderson Bear Claw Retractor

SKU RM-ABCR-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
3 sharp claws, 4 ¾" (12.1 cm)	RM-ABCR-1-1	\$49
5 sharp claws, 4 ¾" (12.1 cm)	RM-ABCR-1-2	\$59

PRODUCT DESCRIPTION

The Anderson Bear Claw Retractor is a specialized handheld plastic surgery instrument engineered to provide exceptional grip and steady traction on cutaneous and subcutaneous tissue flaps. Designed with distinctively curved, sharp claw prongs, this instrument allows surgeons to grasp fibrous tissue planes securely without crushing fragile skin edges.

Optimized for delicate facial rejuvenation and reconstructive procedures, this Anderson Bear Claw Retractor model features a compact overall length of 4 ¾ inches (12.1 cm). It is available in two specialized claw variations: a narrow 3-prong sharp layout (\$49) for targeted digital or periorbital access, and a broader 5-prong sharp layout (\$59) designed to distribute retractive force evenly across larger cutaneous flaps. Its ergonomically contoured handle features a wide thumb-rest to ensure non-slip manual control.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the Anderson Bear Claw Retractor non-sterile and fully autoclavable in protective packaging, delivering high tactile sensitivity and complete corrosion resistance for plastic and aesthetic surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)