



RETRACTOR MAKER

PRODUCT BROCHURE OF

Schuknecht Fisch Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Schuknecht Fisch Retractor

SKU RM-SFR-13-1

\$114 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
2 Teeth (2:2 Array) — 13 cm length, semi-sharp 2:2 prongs, detachable	RM-SFR-13-1-1	\$114
3 Teeth (3:3 Array) — 13 cm length, semi-sharp 3:3 prongs, detachable	RM-SFR-13-1-2	\$124

PRODUCT DESCRIPTION

The Schuknecht Fisch Retractor is an internationally acclaimed otologic and lateral skull base surgical retractor engineered specifically to provide stable, hands-free soft-tissue and periosteal retraction in confined surgical fields. Designed with an articulated self-retaining chassis, this instrument holds postauricular incisions open under calibrated mechanical tension, eliminating assistant hand fatigue while preserving an open visual corridor under operating microscopes.

Built for complete procedural reliability and reprocessing ease, this Schuknecht Fisch Retractor configuration features an overall tool length of 13 cm. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates a semi-sharp prong configuration available in two selectable variations: a 2 teeth (2:2 prong array) layout (\$114) for compact incisions or a 3 teeth (3:3 prong array) layout (\$124) for broader wound margins. Its reprocessing design is built with fully detachable parts, allowing the mechanical chassis to break down quickly for convenient, thorough cleansing and absolute sterilization reliability.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Schuknecht Fisch Retractor to withstand repetitive hospital autoclave cycles, delivering exceptional mechanical stability, corrosion resistance, and complete autoclave durability for otolaryngology suites and neurosurgical operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)