



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hajek Cheek Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hajek Cheek Retractor

SKU RM-HCR-175-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Hajek Cheek Retractor	17.5 cm overall length, atraumatic cheek blade, comfort grip — 1	\$24
Protective Sleeve / Box	Individual transit packaging — 1	\$24

PRODUCT DESCRIPTION

The Hajek Cheek Retractor is an essential oral and maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic displacement of the cheeks, lips, and oral mucosa during dental, alveolar, and orthognathic surgery. Designed with an atraumatic blade profile, this instrument shields soft tissues from surgical drills and rotary burrs while maintaining a wide operative view over the dentition and jaws.

Built to ensure precise operational handling, this Hajek Cheek Retractor configuration features an overall tool length of 17.5 cm (\$24) with a lightweight instrument mass of 0.080 kg that minimizes user fatigue. Its working end incorporates a smooth, atraumatic blade layout that distributes manual counter-traction evenly across buccal tissues, preventing localized pressure ischemia or mucosal tearing. The ergonomic comfort-grip handle provides surgical assistants with steady, non-slip manual control throughout lengthy oral procedures.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Hajek Cheek Retractor to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for dental, oral surgery, and maxillofacial clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)