



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hargis Medial Ramus Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hargis Medial Ramus Retractor

SKU RM-HMRR-1

\$59 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Left-Side Facing — 7 ¾" (19.7 cm)	RM-HMRR-1-1	\$59
Right-Side Facing — 7 ¾" (19.7 cm)	RM-HMRR-1-2	\$59

PRODUCT DESCRIPTION

The Hargis Medial Ramus Retractor is a specialized orthognathic surgical retractor engineered specifically to provide safe, unobstructed access to the medial surface of the mandibular ascending ramus. Designed to slide gently between the medial ramus cortex and the medial pterygoid muscle, this instrument shields the inferior alveolar nerve, artery, and vein above the mandibular foramen during horizontal osteotomy cuts.

Built to accommodate bilateral surgical approaches, the Hargis Medial Ramus Retractor provides a fixed overall length framework of 7 ¾ inches (19.7 cm, \$59). The instrument is manufactured in two mirrored anatomical shapes selectable from the menu layout: a Left-side facing variation and a Right-side facing variation. The curved, atraumatic blade tip wraps snugly around the posterior and medial border of the ramus, preventing tissue slippage while keeping the surgical assistant's hands safely outside the operative line of sight.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Hargis Medial Ramus Retractor to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and facial reconstructive surgeons worldwide with exceptional durability, corrosion resistance, and clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)