



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hoyos Proximal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hoyos Proximal Retractor

SKU RM-HPR-20-1

\$499 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Hoyos Proximal Retractor	20 cm overall length, TC-integrated handle, proximal curved blade	\$499
Protective Sleeve / Box	Individual transit packaging	\$499

PRODUCT DESCRIPTION

The Hoyos Proximal Retractor is an advanced aesthetic plastic surgery retractor engineered specifically to provide superior proximal flap elevation and expansive visual clearance during high-definition body contouring and lipo-abdominoplasty. Designed to handle thick, vascularized tissue planes, this instrument allows surgeons to elevate the epigastric flap toward the xiphoid process with minimal tissue trauma, exposing the upper rectus fascia for precise muscular definition.

Built to optimize field clearance and withstand heavy pulling leverage, this Hoyos Proximal Retractor configuration features an overall tool length of 20 cm (\$499). Its handle assembly integrates a heavy-duty grip framework equipped with a proprietary tungsten carbide handle mechanism that delivers unmatched balance, non-slip manual control, and reduced assistant hand fatigue during multi-hour body sculpting procedures. The polished blade architecture distributes counter-traction evenly to safeguard dermal vascularity.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with tungsten carbide inserts and an anti-glare satin finish, Retractor Maker produces the Hoyos Proximal Retractor to strict clinical standards, providing cosmetic and plastic surgery centers worldwide with exceptional mechanical durability, unyielding precision, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)