



RETRACTOR MAKER

PRODUCT BROCHURE OF

Thompson Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Thompson Retractor System

SKU RM-TRS-1

\$2,499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
2	Elite Rail Clamp	22" (55.9 cm) heavy-duty side rail mount
1	Double Angled Arm	22" (7" × 9" × 7" / 17.8 × 22.9 × 17.8 cm)
1	Lateral Angled Arm	19" (12" × 7" / 30.5 × 17.8 cm) @ 45° (3/8" diameter)
2	Lateral Angled Arms	17" (10" × 7" / 25.4 × 17.8 cm) @ 45° (3/8" diameter)
1	Bilateral Crossbar (Hinged)	34" (11½" × 8½" × 11½" / 29.2 × 21.6 × 29.2 cm)
3	Cam Clip-On	Quick-connect frame sliding clips
1	Micro-Adjustable Handle	10" (25.4 cm) with precise tension screw
3	Balfour Blades with Lips	83 × 73 mm, 83 × 123 mm, 83 × 175 mm
3	Wide Balfour Blades	Interchangeable, broad muscle wall hold
1	Harrington (Sweetheart) Blade	114 × 254 mm (4½" × 10") for deep organ hold
9	Kelly Interchangeable Blades	3" × 3½" (7.6 × 8.9 cm) standard right-angle
1	Malleable Interchangeable Blade Set	25×190 mm, 51×254 mm, 76×254 mm, 102×254 mm
2	St. Marks Interchangeable Blades	51 × 152 mm (2" × 6") with deep lip edge
1	Storage Box	Wooden storage case with customized padding slots

PRODUCT DESCRIPTION

The Thompson Retractor System is a dual-rail, table-mounted self-retaining retractor platform engineered to eliminate assistant fatigue and maintain rock-solid operative visibility in deep anatomical cavities. Built for demanding surgical disciplines, its modular framework anchors securely to both operating table rails using two heavy-duty Elite clamps, supporting an independent bilateral hinged crossbar and multiple angled articulated arms that distribute counter-traction evenly across wide surgical wounds.

Engineered to handle the heavy retraction loads required in liver transplants, open vascular reconstructions, hepatobiliary resections, and radical pelvic dissections, the Thompson Retractor System incorporates quick-connect Cam Clip-On attachments alongside fine-threaded micro-adjustable handles. The system is paired with a diverse surgical blade inventory, including lipped and wide Balfour blades, deep visceral Harrington Sweethearts, St. Marks pelvic blades, standard Kelly retractors, and multi-width malleable ribbons that adapt smoothly to irregular patient anatomy.

Forged exclusively from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a specialized anti-glare satin finish, Retractor Maker supplies the Thompson Retractor System non-sterile and fully autoclavable inside a tailored wooden storage chest designed for organized transport and surgical tray preparation.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)