



RETRACTOR MAKER

PRODUCT BROCHURE OF

Collateral Ligament Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Collateral Ligament Retractors

SKU RM-CLR-1

\$104 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Collateral Ligament Retractor	8" (20.3 cm) overall length, 5 mm & 11 mm blade profiles	\$104
Protective Sleeve / Box	Individual transit packaging	\$104

PRODUCT DESCRIPTION

The Collateral Ligament Retractors are specialized orthopedic surgical instruments engineered specifically to shield and retract the medial collateral ligament (MCL) and lateral collateral ligament (LCL) during total knee replacement. Designed to seat snugly against the tibial plateau and femoral condyle cortices, these instruments keep ligamentous tissues isolated from oscillating bone saws and rotary burrs.

Built for precise joint margin protection, this Collateral Ligament Retractors model features a uniform 8-inch (20.3 cm) overall chassis length (\$104). The working end integrates dual-profile blade widths—incorporating a narrow 5 mm profile for tight joint spaces alongside an 11 mm broad profile for wider ligamentous coverage. Its contoured blade tip hooks securely over the bone edge, preventing accidental slippage during aggressive tibial or femoral resections.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces Collateral Ligament Retractors to withstand heavy mechanical loads, providing orthopedic knee reconstructive surgeons worldwide with exceptional tactile control, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)