



RETRACTOR MAKER

PRODUCT BROCHURE OF

Heaney Simon Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Heaney Simon Vaginal Retractor

SKU RM-HSVR-1

\$139 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
27.0 cm (10.6") — Small, Medium, Large Profiles	RM-HSVR-1-1	\$139
29.2 cm (11.5") — Small, Medium, Large Profiles	RM-HSVR-1-2	\$149

PRODUCT DESCRIPTION

The Heaney Simon Vaginal Retractor is an internationally acclaimed gynecological and pelvic reconstructive surgical retractor engineered specifically to provide deep, stable vaginal wall displacement and apical vault visualization. Designed with an elongated curved spatula blade, this instrument depresses the posterior vaginal wall and pelvic floor smoothly, establishing an expansive visual corridor to the cervix and peritoneal reflection.

Built for heavy pelvic counter-traction, the Heaney Simon Vaginal Retractor is available across two overall lengths: 27 cm (10.6 inches, \$139) and an extended 29.2 cm (11.5 inches, \$149) reach for deep pelvic or parous patient profiles. The instrument is categorized across three standard operational size profiles via the configuration menu, ensuring optimal blade surface area for variable vaginal capacities. Its ergonomic drop handle provides surgical assistants with steady pulling leverage, minimizing hand strain during prolonged procedures.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Heaney Simon Vaginal Retractor to withstand repetitive hospital autoclave cycles, providing obstetrical, gynecological, and urogynecological teams worldwide with exceptional durability, corrosion resistance, and clinical reliability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)