



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Markham-Meyerding Hemilaminectomy Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Markham-Meyerding Hemilaminectomy Retractor

**SKU** RM-MMHR-17-1

**\$124 USD**

## AVAILABLE SIZES & PRICING

| SIZE / OPTION                | SKU / DETAILS  | PRICE (USD) |
|------------------------------|----------------|-------------|
| Small — 40 mm x 25 mm blade  | RM-MMHR-17-1-1 | \$104       |
| Medium — 50 mm x 25 mm blade | RM-MMHR-17-1-2 | \$114       |
| Large — 60 mm x 25 mm blade  | RM-MMHR-17-1-3 | \$124       |

## PRODUCT DESCRIPTION

The Markham-Meyerding Hemilaminectomy Retractor is a specialized spinal surgical retractor engineered specifically to establish and maintain an unyielding unilateral operative window during posterior lumbar hemilaminectomy and discectomy. Designed with an asymmetrical working framework, this retractor anchors against the spinous process while simultaneously retracting heavy paraspinal musculature laterally away from the laminar corridor.

Built to clinical framework specifications, this Markham-Meyerding Hemilaminectomy Retractor retains a fixed overall length of 17 cm (\$124) while offering two mirrored anatomical layouts selectable via the purchase configuration menu. Option 1 features a 63.5 mm long x 25 mm wide blade on the left side paired with a single hook on the right side, while Option 2 places the blade on the right side paired with a left hook. Its multi-position mechanical ratchet lock sustains continuous counter-traction automatically, eliminating assistant hand fatigue throughout deep canal decompression.

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Markham-Meyerding Hemilaminectomy Retractor to withstand repetitive hospital autoclave cycles, providing spine centers, neurosurgical suites, and orthopedic trauma units worldwide with unyielding structural rigidity, high corrosion resistance, and complete clinical durability.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)