



RETRACTOR MAKER

PRODUCT BROCHURE OF

Murphy Ball Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Murphy Ball Retractor

SKU RM-MBR-35-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
2.0 mm Spherical Tip	RM-MBR-35-1-1	\$39
3.0 mm Spherical Tip	RM-MBR-35-1-2	\$39
4.0 mm Spherical Tip	RM-MBR-35-1-3	\$39

PRODUCT DESCRIPTION

The Murphy Ball Retractor is an indispensable handheld surgical instrument engineered specifically to provide safe, non-traumatic tissue exploration, nerve manipulation, and vessel displacement in deep surgical cavities. Designed with an ultra-smooth spherical ball tip, this retractor hooks gently around delicate tubular structures, allowing surgeons to mobilize tissue planes without risking accidental puncture, tearing, or intimal injury.

Built for extended deep-cavity reach, this Murphy Ball Retractor maintains a uniform overall instrument length of 35 cm (\$39) that provides exceptional clearance in deep pelvic, spinal, and abdominal corridors. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry features a 60-degree angled hook design paired with selectable ball tip diameters: 2.0 mm, 3.0 mm, or 4.0 mm rounded spherical tips designed to facilitate gentle tissue manipulation while minimizing potential puncture trauma. Its balanced handle delivers refined tactile feedback.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Murphy Ball Retractor to withstand repetitive hospital autoclave cycles, providing general, spinal, and vascular operating theaters worldwide with exceptional corrosion resistance, unyielding structural rigidity, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)