



RETRACTOR MAKER

PRODUCT BROCHURE OF

Farabeuf Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Farabeuf Retractor

SKU RM-FBR-SET-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
12 cm (approx. 4.7") set	RM-FBR-SET-1-1	\$24
15 cm (approx. 6.0") set	RM-FBR-SET-1-2	\$29

PRODUCT DESCRIPTION

The Farabeuf Retractor is an internationally recognized double-ended surgical retractor engineered specifically to provide reliable, atraumatic wound margin and soft-tissue retraction during general, orthopedic, and plastic surgical procedures. Designed with contrasting right-angled blade lengths on opposite ends, this instrument allows surgical teams to transition effortlessly between shallow subcutaneous tissue retraction and deeper fascial elevation.

Built for clinical versatility across varying surgical depths, this Farabeuf Retractor set is available in two primary length options: 12 cm (\$24) for shallow, localized incisions and 15 cm (\$29) for deeper operative fields. Each retractor features a short right-angled blade on one end and a slightly longer right-angled blade on the reverse end, both finished with rounded, non-traumatic margins that prevent tissue bruising. The solid flat central handle provides surgical assistants with steady, non-slip manual control.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Farabeuf Retractor to withstand repetitive hospital sterilization cycles, providing operating theaters worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)